

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sephadex® G-50 - Fine

Product Number: **G5080**
CAS Number: 9048-71-9
MDL: MFCD00132771

TEST	Specification
------	---------------

Loss on Drying	≤ 10 %
----------------	--------

Particle Size Distribution	≥ 80 %
----------------------------	--------

Within 20 - 80 microns

Note

Sephadex is a registered trademark of Cytiva

Specification: PRD.3.ZQ5.10000015085

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

