

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Silicone oil - Dow Corning Corporation 200® fluid, viscosity 60,000 cSt (25 °C)

Product Number:

181838


CAS Number:

63148-62-9

MDL:

MFCD00132673

Formula:

(-Si(CH<sub>3</sub>)<sub>2</sub>O-)n**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Viscous Liquid

Infrared spectrum

Conforms to Structure

Loss on Drying

≤ 2 %

Viscosity

57000 - 63000 cps

Flash Point

≥ 610 °F

Product of Supplier

Confirmed

Dow Corning

Specification: PRD.1.ZQ5.10000029356

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.