

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

(–)-1,2-Bis[(2R,5R)-2,5-diethylphospholano]benzene – kanata purity

Product Number:

668494

CAS Number:

136705-64-1

MDL:


MFCD00142335

Formula:

C<sub>22</sub>H<sub>36</sub>P<sub>2</sub>

Formula Weight:

362.47 g/mol



### TEST

### Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Proton NMR spectrum

Conforms to Structure

Phosphorus NMR Spectra

Conforms to Structure

Carbon

68.5 - 77.3 %

Optical Rotation

-220.0 - -200.0 deg

c = 1%; Chloroform

Specification: PRD.1.ZQ5.10000044571

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.