

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Cobalt - powder, <150 µm, ≥99.9% trace metals basis

Product Number:**266647**

CAS Number:

7440-48-4

MDL:

MFCD00010935

Formula:

Co

Formula Weight:

58.93 g/mol

Co

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Grey to Dark Brown

Appearance (Form)

Powder

Size

Confirmed

-100 Mesh

ICP Major Analysis

Confirmed

Confirms Cobalt Component

Purity

Conforms

> = 99.9% Based On Trace Metals Analysis

Trace Metal Analysis

≤ 1000.0 ppm

Specification: PRD.0.ZQ5.10000088406

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.