

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Glass beads, acid-washed – 425–600 µm (30–40 U.S. sieve)

Product Number: G8772

TEST	Specification
Appearance (Form) May contain dark particles/beads	Beads
Appearance (Color) Gray, blue, green or brown	Conforms
Size 425-600 microns	Conforms

Specification: PRD.2.ZQ5.10000010103

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.