

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
ProteoMass™ Angiotensin II MALDI-MS Standard - vial of 10 nmol

Product Number: **A8846**
CAS Number: 68521-88-0
MDL: MFCD00161484
Formula: C₅₀H₇₁N₁₃O₁₂
Formula Weight: 1,046.18 g/mol

H-Asp-Arg-Val-Tyr-Ile-His-Pro-Phe-OH

TEST	Specification
------	---------------

Suitability	Suitable
Suitable for use as a standard for MALDI Mass Spectrometry.	
ug protein/vial	8.4 - 18.8
by Micro BCA Protein Assay	
nmoles/vial	8.0 - 18.0
by Micro BCA Protein Assay	

Specification: PRD.0.ZQ5.10000012444

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.