


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(3-Bromopropyl)trimethylammonium bromide - 97%

Product Number: 347604
CAS Number: 3779-42-8
MDL: MFCD00075396
Formula: C₆H₁₅Br₂N
Formula Weight: 261.00 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Yellow
Appearance (Form)	Solid or Powder or Crystals or Chunks
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃ After Oxygen Combustion	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000055287

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.