

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
SODIUM CARBONATE, ANHYDROUS, FREE-FLOWING

Product Number: **900505**

CAS Number: 497-19-8

Formula: Na_2CO_3

Formula Weight: 105.99 g/mol

TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, H ₂ O	
Titration with HCl	≥ 99.5 %
Dried Basis	
Recommended Retest Period	-----
5 Years	

Specification: PRD.0.ZQ5.10000073388

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.