

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Cytochalasin D – from *Zygosporium mansonii*, ≥98% (TLC and HPLC), powder

Product Number: **C8273**  
CAS Number: 22144-77-0

Formula: C<sub>30</sub>H<sub>37</sub>NO<sub>6</sub>  
Formula Weight: 507.62 g/mol  
Storage Temperature: -20 °C



TEST	Specification
------	---------------

Appearance (Form)	Powder
Appearance (Colour)	White to White w/ Yellow Cast
Solubility (Solvent)	Chloroform
Solubility (Conc)	1.90 - 2.10 mg/ml
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless to Faint Yellow
Purity (TLC)	≥ 98.00 %
Elemental Anal. (%C anhydrous)	70.48 - 71.48 %
Elemental Anal. (%N anhydrous)	2.46 - 3.06 %
NMR (Solvent)	DMSO-d <sub>6</sub>
Identity by NMR	Consistent
Purity (HPLC)	≥ 98.00 %
Recommended Retest Period	-----
3 Years	

Specification: PRE.1.ZQ5.10000001403

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.