

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Collagen from calf skin - Bornstein and Traub Type I (Sigma Type III), solid

Product Number:

C3511

CAS Number:

9007-34-5

MDL:

MFCD00130825

Storage Temperature:

2 - 8 °C

TEST**Specification**

Appearance (Color)

White

Appearance (Form)

Solid

Solubility

Conforms

Opalescent Viscous Solution at 5 mg/mL in H₂O, pH 3.0 with acetic acid

Water (by Karl Fischer)

≤ 15 %

Suitability

Conforms

Suitable for use as a substrate for Collagenase

Recommended Retest Period

8 years

Specification: PRD.0.ZQ5.10000016072

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.