


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Naphthol AS-MX phosphate

Product Number: 855
CAS Number: 1596-56-1
MDL: MFCD00042711
Formula: C₁₉H₁₈NO₅P
Formula Weight: 371.32 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Solution
Appearance (Turbidity)	Clear
Essentially free of precipitate	
Expiration Date	Conforms
30 months	
Performance	Conforms
When staining with Sigma-Aldrich Procedure No. 85, cells stain in accordance with the package insert.	

Specification: PRD.1.ZQ5.10000038276

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.