


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Tetrabutylammonium fluoride solution - 75 wt. % in H<sub>2</sub>O

Product Number: **361399**



MDL: MFCD00011747

Formula: C<sub>16</sub>H<sub>36</sub>FN

Formula Weight: 261.46 g/mol

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Light Yellow to Light Brown
Appearance (Form)	Viscous Liquid
Proton NMR spectra	Conforms to Structure
Purity (Titration by HCL04)	70 - 80 %
Water (by Karl Fischer)	20 - 30 %

Specification: PRD.0.ZQ5.10000068984

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.