

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:Protein Kinase C β 1 isozyme human - \geq 95% (SDS-PAGE), recombinant, expressed in baculovirus infected insect cells, buffered aqueous glycerol solution**Product Number:****P1787****MDL:**

MFCD03458220

Storage Temperature:

-70 °C

TEST	Specification
mg protein/ml (Bradford)	\geq 0.1
units/mg protein	\geq 700
One unit will transfer 1 nmol of phosphate to histone H1 in 1 min at pH 7.4 at 30 °C	
Purity by SDS Electrophoresis	\geq 95 %

Specification: PRD.0.ZQ5.10000015515

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.