

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Peroxidase from horseradish - Type XII, essentially salt-free, lyophilized powder, ≥ 250 units/mg solid (using pyrogallol)

Product Number:

P8415

CAS Number:

9003-99-0

MDL:

MFCD00071339

Storage Temperature:

2 - 8 °C

| TEST | Specification |
|------|---------------|
|------|---------------|

units/mg solid

 ≥ 250 Unit definition: One unit will form 1.0 mg
Purpurogallin from Pyrogallol in 20 seconds
at pH 6.0 at 20 deg C.

RZ

 ≥ 3.0

Ratio

Specification: PRD.0.ZQ5.10000069399

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.