

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Albumin from chicken egg white - lyophilized powder, ≥90% (agarose gel electrophoresis)

 Product Number:
 A5378

 CAS Number:
 9006-59-1

 MDL:
 MFCD00130424

Storage Temperature: 2 - 8 ℃

TEST Specification

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)
50 mg/ml in H2O
Water (by Karl Fischer)
Nitrogen Content

Agarose Electrophoresis

Albumin

Recommended Retest Period

3 years

Specification: PRD.1.ZQ5.10000064103

White to Yellow

Pow der

Faint Yellow to Yellow Clear to Slightly Hazy

≤ 5 %

13.0 - 16.0 %

> 90 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.