

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hexamethonium chloride - solid (hygroscopic)

Product Number:**H2138**

CAS Number:

60-25-3

MDL:

MFCD00031563

Formula:

C₁₂H₃₀Cl₂N₂

Formula Weight:

273.29 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, H ₂ O	
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 13 %
Purity (TLC)	≥ 99 %
Purity Titr. by HCL04 anhyd.	≥ 98 %

Specification Date : 10/29/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.