


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cyclopropylboronic acid

Product Number: 597988
CAS Number: 411235-57-9
MDL: MFCD04038750
Formula: C₃H₇BO₂
Formula Weight: 85.90 g/mol
Storage Temperature: -20 °C



TEST	Specification
------	---------------

Appearance (Color)	White to Light Yellow to Light Orange and White to Light Grey
Appearance (Form)	Powder
Proton NMR Spectrum	Conforms to Structure
Boron NMR Spectrum	Conforms to Structure
Purity (Titration by NaOH)	≥ 90.0 %
Note	Note
May Contain 3-5% Cyclopropanol	
May Contain 5-10% Boric Acid	

Specification: PRD.3.ZQ5.10000035553

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.