

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-phospho-Synapsin I (pSer⁶⁰³) antibody produced in rabbit - ~1 mg/mL, affinity isolated antibody, buffered aqueous solutionProduct Number: **S8192**

MDL: MFCD03792019

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
Solution in 10 mM HEPES, pH 7.5, containing 150 mM NaCl, 100 ug/mL BSA, and 50% Glycerol	
Working Dilution	Conforms
by Direct Immunoblotting	
1:1,000 using rat brain hippocampal tissue homogenate	

Specification: PRD.0.ZQ5.10000040091

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.