

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Myosin IIB, non-muscle antibody produced in rabbit - affinity isolated antibody, buffered aqueous solution

Product Number:

M7939

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	1.0 - 1.5 mg/ml
WB-Cell Line/Tissue Extract	MDCK
WB Titer	1:200-1:400
WB-Cell Line/Tissue Extract	Jurkat Cells
WB-Western Blot Titer	1:200-1:400
WB-Cell Line/Tissue Extract	NRK rat kidney cells
WB-Western Blot Titer	1:200-1:400
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000019780

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.