


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tetraiodomethane - 97%

Product Number: 269220
CAS Number: 507-25-5
MDL: MFCD00001067
Formula: CI₄
Formula Weight: 519.63 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Conforms to Requirements
Red to Purple	
Appearance (Form)	Powder or Solid or Crystals
X-Ray Diffraction	Conforms to Structure
Titration by AgNO ₃	96.5 - 103.5 %
After Oxygen Combustion	

Specification: PRD.0.ZQ5.10000036087

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.