

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Human Lambda Light Chains (Bound and Free) antibody produced in mouse - clone HP-6054, ascites fluid

Product Number: **L6522**

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Faint Brown to Very Dark Brown to Very Dark Red
Protein by Biuret	30 - 60 mg/ml
Mouse IgG2a by RID	4.0 - 7.0 mg/ml
ELISA - Titer	1:1000-1:2000
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000019441

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.