

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Caffeine-(trimethyl-d₉) - 99 atom % D, 98% (CP)

Product Number:

725625

CAS Number:


72238-85-8

Formula:

C₈D₉H_N4O₂

Formula Weight:

203.16 g/mol



| TEST | Specification |
|-----------------------|---------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder |
| Purity (HPLC) | ≥ 98 % |
| Initial Melting Point | |
| Final Melting Point | |
| D Enrichment | ≥ 99 atom % |

Specification: PRD.0.ZQ5.10000102900

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

