

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification


Product Name: (S)-N-Methyl-N-diphenylphosphino-1-[(R)-2-(diphenylphosphino)ferrocenyl]ethylamine

Product Number: **682314**

CAS Number: 406681-09-2

Formula: C₃₇H₃₅FeNP₂

Formula Weight: 611.5 AMU



TEST	Specification
Appearance (Color)	Yellow to Orange
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Carbon	69.8- 75.6 %
Optical Rotation c = 1%; Chloroform	273.0- 283.0 deg

Remarks:

Specification Date : 12/17/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.