

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium chloride – tablet

Product Number: **S6150**
CAS Number: 7647-14-5
MDL: MFCD00003477
Formula: NaCl
Formula Weight: 58.44 g/mol

NaCl

TEST	Specification
------	---------------

Total Tablet Weight	1615 - 1785 mg
Appearance	White
Diameter	5/8 inches
Form	Round Tablet
Dissolution Time	≤ 10 min
1 tablet in 100 mL purified water	
pH	5.5 - 6.5
1 tablet in 100 mL purified water	

Specification: PRD.2.ZQ5.10000006757

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.