

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

ProteoMass™ Aldolase MALDI-MS Standard – vial of 10 nmol, average mol wt 39,211.28 Da by calculation

Product Number:

A9096

CAS Number:

9024-52-6

MDL:

MFCD00130453

TEST	Specification
------	---------------

Protein (nanomoles/vial)

8.0 - 18.0

by UV Absorbance

ug protein/vial

313.7 - 705.8

by UV Absorbance

Suitability

Suitable

Suitable for use as a standard for MALDI Mass Spectrometry

Specification: PRD.0.ZQ5.10000066208

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.