

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Chemiluminescent Peroxidase Substrate-1 – sufficient for 60 mL substrate

Product Number: **CPS160**

MDL: MFCD03456389

TEST	Specification
Kit Comp (Component 1) C9232	-----
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Solution
pH	7.25 - 7.75
Kit Comp (Component 2): C9107	-----
Appearance (Turbidity)	Clear
Appearance (Color)	Conforms
Colorless to pinkish tint	
Appearance of Solution (Form)	Liquid
pH	9.35 - 9.85

Specification: PRD.1.ZQ5.10000009361

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.