

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Leukemia Inhibitory Factor antibody produced in goat – affinity isolated antibody, lyophilized powder

Product Number:

L1044

Storage Temperature:

-20 °C

### TEST

### Specification

Appearance (Form)

from a 0.2 µm filtered solution in PBS with 5% trehalose

Endotoxin Assay

Western Blot

Neutralization of Bioactivity

10-60 ng/ml

Lyophilized Powder

Less than 0.1 Eu per 1 µg

0.1 to 0.2 µg/ml

Conforms

Specification: PRD.2.ZQ5.10000011584

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

