


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-Vinyl-2-pyrrolidinone - contains sodium hydroxide as inhibitor, ≥99%

Product Number: V3409
CAS Number: 88-12-0
MDL: MFCD00003197
Formula: C₆H₉NO
Formula Weight: 111.14 g/mol



TEST

Specification

| | |
|--|---------------------------|
| Appearance (Color) | Colorless to Light Yellow |
| Material may darken in color over time | |
| Appearance (Form) | Liquid |
| Refractive index at 20 °C | 1.510 - 1.514 |
| Infrared spectrum | Conforms to Structure |
| Purity (GC) | ≥ 99.0 % |
| Inhibitor | Present |
| Confirms Na Component | |

Specification: PRD.2.ZQ5.10000024064

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.