

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

**Product Specification**

Product Name:

9-Fluorenyl methacrylate 97%

Product Number:

**733342**

CAS Number:

46969-53-3

Formula:


C<sub>17</sub>H<sub>14</sub>O<sub>2</sub>

Formula Weight:

250.29 g/mol

Storage Temperature:

2 - 8 °C



**TEST**

**Specification**

Appearance (Color)	White to Tan
Appearance (Form) Powder or Crystalline Powder	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	78.7- 84.5 %
Purity (GC)	≥96.5 %
Wavelength Report Results	
Extinction Coefficient Report Results	

**Remarks:**

Specification Date : 12/08/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.