

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**Product Name:
N-Boc-tyramine 97%**Product Number:****691550**

CAS Number:


64318-28-1

Formula:

C₁₃H₁₉NO₃

Formula Weight:

237.29 g/mol



TEST	Specification
Appearance (Color)	White to Tan
Appearance (Form)	Conforms to Requirements
Solid, Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	63.5- 68.1 %
Nitrogen	5.7- 6.1 %
Purity (HPLC)	≥96.5 %

Remarks:

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.