

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Salinomycin, Ready Made Solution - from *Streptomyces albus*, ≥98% (HPLC)

Product Number:

S6201



Formula:

C₄₂H₇₀O₁₁

Formula Weight:

751.00 g/mol

Storage Temperature:

-20 °C

| TEST | Specification |
|------|---------------|
|------|---------------|

Appearance (Form)

Clear Solution

Appearance (Colour)

Colorless

Concentration

2.00 - 2.00 mg/ml

Purity (HPLC)

≥ 98.00 %

Recommended Retest Period

2 Years

Specification: PRE.0.ZQ5.10000004407

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.