

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Gelatin from porcine skin - gel strength 300, Type A

Product Number: G2500
CAS Number: 9000-70-8
MDL: MFCD00081638

TEST	Specification
Appearance (Color)	Off-White to Yellow
Appearance (Form)	Powder
Solubility (Color)	Faint Yellow
Solubility (Turbidity)	Clear to Hazy
50 mg/ml, H2O	
% Protein (Biuret)	70 - 90
Gel Strength (Bloom No.)	275 - 325
(approx 300)	
Recommended Retest Period	-----
5 years	

Specification: PRD.0.ZQ5.10000029179

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.