


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Imino-1-imidazolidineacetic acid - 98%

Product Number: 377627
CAS Number: 35404-50-3
MDL: MFCD00134455
Formula: C₅H₉N₃O₂
Formula Weight: 143.14 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	40.9 - 43.0 %
Nitrogen	28.6 - 30.1 %
Purity (TLC)	≥ 98 %
Solubility (Turbidity) 25 mg/ml, 1N HCl	Clear
Solubility (Color)	Colorless

Specification: PRD.0.ZQ5.10000066397

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.