# SIGMA-ALDRICH®

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

# **Product Specification**

Product Name: Spectinomycin dihydrochloride pentahydrate - powder, BioReagent, suitable for cell culture

#### **Product Number:**

CAS Number: MDL: Formula: Formula Weight: Storage Temperature: **S4014** 22189-32-8 MFCD00150886 C14H24N2O7 · 2HCI · 5H2O 495.35 g/mol 2 - 8 ℃

# 

## TEST

Appearance (Color) Appearance (Form) Solubility (Color) Solubility (Turbidity) 50 mg/ml, H2O Infrared spectrum Water (by Karl Fischer) µg/mg Cell Culture Test Endotoxin Level Titration with AgNO3 Recommended Retest Period 5 years

Specification: PRD.0.ZQ5.10000042337

### Specification

White to Off-White Pow der Colorless to Light Yellow Clear Conforms to Structure 16.5 - 19.5 % ≥ 603 Pass ≤ 10 EU/mg 98.0 - 102.0 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.