

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(tetrahydrofuran) - average $M_n \sim 250$

Product Number:

345261


CAS Number:

25190-06-1

MDL:

MFCD00148879

Formula:

 $H(OCH_2CH_2CH_2CH_2)_nOH$ **TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Viscous Liquid

Infrared spectrum

Conforms to Structure

Molecular Number

230 - 270

Color Test

 ≤ 40 APHA

Hydroxyl Value, mg KOH/g

415.6 - 487.8

Specification Date : 06/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.