

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

StratoSpheres™ PL-NCO (Isocyanate) resin - 50-100 mesh, extent of labeling: 1.5-2.0 mmol/g loading, 1% cross-linked

Product Number:**547530****CAS Number:**


59990-69-1

MDL:

MFCD01867272

Formula:C₂₆H₂₅NO**Formula Weight:**

367.48 g/mol

**TEST****Specification**

Appearance (Color)	Tan
Appearance (Form)	Beads
Infrared spectrum	Conforms to Structure
Particle Size	Confirmed
150 - 300 Microns	
Particle Size (Median)	150 - 300 micron
Particle Size	
Volume Between 150-300 Microns	
Swelling Volume	
Solvent: Dichloromethane	
Loading, Molecular	≥ 1.5 mmol/g
Residual Solvent	Conforms
Dried to Constant Weight	
Note	Conforms
StratoSpheres is a trademark of Agilent Technologies	

Specification: PRD.0.ZQ5.10000012405

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.