

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,6-Diphenyl-1,3,5-hexatriene - 98%

 Product Number:
 D208000

 CAS Number:
 1720-32-7

 MDL:
 MFCD00004793

Formula: C18H16 Formula Weight: 232.32 g/mol

TEST Specification

Appearance (Color)
Appearance (Form)

Powder, Flakes or Crystals

Infrared spectrum Purity (HPLC)

Solubility (Turbidity) c = 10 mg/ml, CHCI3

Solubility (Color)

Yellow to Tan

Conforms to Requirements

Conforms to Structure

≥ 97.5 %

Clear to Hazy

Yellow to Yellow-Green

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.