

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Rhodium(III) nitrate solution - ~10 % (w/w) (Rh in >5 wt. % HNOs wt. % HNO<SUB>3

Product Number: 309206

Rh

MDL: MFCD00011201 Formula: N3O9Rh

Formula Weight: 288.92 g/mol

TEST Specification

Appearance (Color) Conforms to Requirements

Yellow to Orange to Brown

Appearance (Form)

Purity (Titration by NaOH)

Liquid

> 5.0 %

wt% HNO3

Rhodium (Rh) 8 - 11 %

Specification: PRD.0.ZQ5.10000038516

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.