

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
EX-CELL® Advanced™ CHO Feed 1 (without glucose)

Product Number: 24368C

Storage Temperature: 2 - 8 °C

TEST	Specification
Origin	Non-Animal Source
Osmolality w/o NaHCO ₃	450 - 550 mOs/kg
Solubility soluble at 1x	Pass
Endotoxin Level	≤ 10 EU/ml
Site of Manufacturing	Country of Origin
Appearance (Color)	Off-White to Tan
Appearance (Form)	Free Flowing Powder

Specification: PRD.0.ZQ5.10000078432

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.