

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Aluminum potassium sulfate dodecahydrate - ACS reagent, ≥98%

Product Number:

237086

AlK(SO₄)₂ · 12H₂O

CAS Number:

7784-24-9

MDL:

MFCD00149143

Formula:

AlK₂O₈S₂ · 12H₂O

Formula Weight:

474.39 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Compleximetric EDTA	98.0 - 102.0 %
X-Ray Diffraction	Conforms to Structure
Insoluble matter	≤ 0.005 %
Chloride (CL)	≤ 5 ppm
Iron (Fe)	≤ 0.001 %
Heavy Metals	≤ 0.001 %
by ICP-OES	
Sodium (Na)	≤ 0.02 %
Ammonium (NH ₄)	≤ 0.005 %
Meets ACS Requirements	Current ACS Specification

Specification: PRD.2.ZQ5.10000058368

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.