


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2'-Bromo-2,6-dimethoxybiphenyl - 97%

Product Number: 655481
CAS Number: 755017-61-9
MDL: MFCD07369791
Formula: C₁₄H₁₃BrO₂
Formula Weight: 293.16 g/mol



TEST

Specification

Appearance (Color)	White to Light Yellow to Light Orange
Appearance (Form)	Powder or Crystals or Powder with Chunk(s) or Crystals with Chunk(s)
Infrared Spectrum	Conforms to Structure
Purity (GC)	≥ 96.5 %
Titration by AgNO ₃ After Oxygen Combustion	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000081815

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.