


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Barium hydroxide solution – 0.3 N

Product Number: **B4059**



MDL: MFCD00003443

Formula: H_2BaO_2

Formula Weight: 171.34 g/mol

TEST	Specification
------	---------------

Appearance	Pass
Clear colorless liquid, may contain some small insoluble particles.	
Normality	0.298 - 0.305 N
by HCl Titration	

Specification: PRD.2.ZQ5.10000009330

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.