

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

L-15 Medium (Leibovitz) - With L-glutamine, liquid, sterile-filtered, suitable for cell culture

Product Number: **L1518**

MDL: MFCD00217482

Storage Temperature: 2 - 8 °C

TEST	Specification
Osmolality	305 - 338 mOs/kg
pH	7.3 - 7.9
Appearance (Turbidity)	Clear
Appearance (Form)	Solution
Sterility by USP guidelines	Pass
Endotoxin Level	≤ 1.0 EU/ml
Galactose Concentration	0.81 - 0.99 g/l
Cell Culture Test	Pass
Cell Line	Attached Cell Line

Specification: PRD.0.ZQ5.10000038317

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.