

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Amberlyst® 15 hydrogen form - wet, moisture ~48%

Product Number:

216399

CAS Number:

39389-20-3

MDL:

MFCD00145841

**TEST****Specification**

Appearance (Color)

Brown to Very Dark Brown and Grey and

Appearance (Form)

Brown-Grey and Grey-Brown

Exchange Capacity

Granules or Pellets or Beads

Registered Trademark

≥ 1.7 meq/mL

Amberlyst is a registered trademark of Rohm and Haas Co.

Confirmed

Specification Date : 04/07/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.