


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(±)-CPP - solid

Product Number: **C104**
CAS Number: 100828-16-8
MDL: MFCD00055136
Formula: C₈H₁₇N₂O₅P
Formula Weight: 252.2 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility	Pass
Soluble in Water	
Proton NMR spectrum	Conforms to Structure
Water Content	≤ 1.4
Expressed in Mole	
Carbon (anhydrous)	37.1 - 39.1 %
Nitrogen (anhydrous)	10.8 - 11.4 %
Purity (TLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000012578

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.