

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Diethylamine-d₁₁ – 98 atom % D, 98% (CP)

Product Number:

617075

CAS Number:


1173019-51-6

Formula:

C₄D₁₁N

Formula Weight:

83.98 g/mol



TEST	Specification
------	---------------

Appearance (Clarity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Purity (GC)

≥ 98 %

Proton NMR Spectrum

Conforms

Proton NMR

≥ 98 %

(ATOM% D)

Specification: PRD.0.ZQ5.10000103159

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

