

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 α -Cyano-4-hydroxycinnamic acid – suitable for MALDI-TOF MS

Product Number:


C8982

CAS Number:

28166-41-8

MDL:

MFCD00004204



Formula Weight:

189.17 g/mol

TEST

Specification

Appearance (Color)

White to Yellow

Form

Powder or Crystals

Suitability

Suitable

THIS LOT TESTED AND FOUND SUITABLE FOR USE IN
MALDI-TOF MASS SPECTROSCOPY

% Purity

≥ 99.0 %

Specification: PRD.1.ZQ5.10000077383

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.