


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hydrocortisone - γ -irradiated, powder, BioXtra, suitable for cell culture

Product Number: H0135
CAS Number: 50-23-7
MDL: MFCD00011654
Formula: C₂₁H₃₀O₅
Formula Weight: 362.46 g/mol



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Turbidity)	Clear
1 mg/mL, EtOH	
Sterility by USP guidelines	Pass
Purity (HPLC)	≥ 80 %
Cell Culture Test	Pass
Sterilized (Gamma Irradiated)	Conforms
Irradiation Dosage	Conforms
25.0 to 35.0 kGy	

Specification: PRD.0.ZQ5.10000036199

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.