

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anti-Mouse IgG (whole molecule) antibody produced in rabbit – IgG fraction of antiserum, buffered aqueous solution

Product Number: M7023

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	5.0 - 12.0 mg/ml
Antibody Concentration-Method	Quantitative Precipitin Assay
Antibody Concentration	2.0 - 4.5 mg/ml
ELISA - Titer	1:80,000-1:120,000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000015909

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.