


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(R)-[(RuCl(T-BINAP))₂(μ-Cl)₃][NH₂Me₂]

Product Number: 692441
CAS Number: 749935-02-2
MDL: MFCD09753034
Formula: C₉₈H₈₈Cl₅NP₄Ru₂
Formula Weight: 1,783.05 AMU



TEST

Specification

Appearance (Color)	Orange to Brown
Appearance (Form)	Powder
Proton NMR Spectrum	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Carbon	62.7 - 69.3 %
Optical Rotation	1400.0 - 1700.0 deg
c = 0.5%; Chloroform	
Product of Supplier	Confirmed
Takasago International Corporation	

Specification: PRD.1.ZQ5.1000004684

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.